

Nano-Check™ SIRS Detection Test PCT / CRP



Nano-Check™ PCT
One-Step quantitative assay in conjunction with Nano-Checker™ 710V Reader

Nano-Check™ CRP
Two-Step quantitative assay in conjunction with Nano-Checker™ 710V Reader

Product Specification

	Catalog No.	Specimen	Storage Condition	Shelf-Life	Turn Around Time	Measuring Range
Nano-Check™ PCT	ND-CD7155	80 µL of whole blood / plasma / serum	RT (2~30°C)	15 months	15 minutes	0.1 ng/mL – 21.00 ng/mL
Nano-Check™ CRP	ND-CD-7151	5 µL of whole blood 10 µL of plasma / serum	RT (2~30°C)	15 months	15 minutes	2 µg/mL – 100 µg/mL